

Raf1 (phospho Ser621) antibody

Cat. No. GTX24767

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, FCM, IP
Reactivity	Human, Rat

References (2)
 Package
 50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	1:20
IP	Assay dependent

Not tested in other applications.

Calculated MW 73 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA, 50% Glycerol
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against a chemically synthesized phosphopeptide derived from a region of human c-Raf that contains serine 621. The sequence is conserved in mouse and rat.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX24767 WB Image

WB analysis of immunoprecipitates prepared HEK293 cells overexpressing wild-type c-Raf and stimulated with EGF (Lane 1-4) or c-Raf mutant S621A (Lane 5) using GTX24767 Raf1 (phospho Ser621) antibody.



For full product information, images and publications, please visit our [website](#).