

Lamin A + C antibody [JoL3]

Cat. No. GTX24789

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
1 ml

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1/50

For ICC/IF: Use at 1:10

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	74 kDa. (Note)
Product Note	Common tail domain of Lamin A & C from human, mouse & rat

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant bacterially expressed human lamin A
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).