

MEK3 antibody

Cat. No. GTX25428

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

For IP and WB: Use at a dilution of 1:1,500 - 1:2,500. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW 39 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein - full length human MEK3.
Purification	Protein A affinity purified From polyclonal serum
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX25428 WB Image

Western blot detection of MKK3 protein in RAW cell (1) and RBL transfected with NIK(2).



For full product information, images and publications, please visit our [website](https://www.genetex.com).