

## A431 (human epidermoid carcinoma) cell line slides

Cat. No. GTX25525

Applications	ICC/IF
Species	Human

References ( 1 )  
Package  
5 slide

## PRODUCT

A-431 (epidermis; epidermoid carcinoma)

## Summary

Human A-431 cells were cultured in DMEM medium and harvested at the log phase of growth. The A-431 cells were fixed in 4% paraformaldehyde and arrayed on a 12-well (5 mm) adhesive coated slide, with each well's surface specifically treated to enhance cellular attachment and to minimize background staining. Each well was loaded with approximately  $7.5 \times 10^3$  A-431 cells. A follow-up heat dehydration process was done to ensure the attachment of A-431 cells and the stability of cellular proteins.

Tissue type	Tumor
-------------	-------

## Applications

## Application Note

A-431 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

## Properties

Storage	Stable for up to 6 months at 4°C under desiccate conditions.
---------	--------------------------------------------------------------

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).