

# EL-4 (mouse T lymphocyte; lymphoma) cell line slides

**Cat. No. GTX25544**

<b>Application</b>	ICC/IF
<b>Species</b>	Mouse

**Package**  
5 slide

## PRODUCT

EL4 (T lymphocyte; lymphoma)

### Summary

Mouse EL4 cells were cultured in DMEM and harvested at the log phase of growth. The EL4 cells were fixed in 4% paraformaldehyde and arrayed on a 12-well (5 mm) adhesive coated slide, with each well's surface specifically treated to enhance cellular attachment and to minimize background staining. Each well was loaded with approximately  $7.5 \times 10^3$  EL4 cells. A follow-up heat dehydration process was done to ensure the attachment of EL4 cells and the stability of cellular proteins.

## APPLICATION

### Application Note

EL4 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

## PROPERTIES

**Storage** Stable for up to 6 months at 4°C under desiccate conditions.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).