

P815 (mouse mast cell; mastocytoma) cell line slides

Cat. No. GTX25547

| | |
|---------------------|--------|
| Applications | ICC/IF |
| Species | Mouse |

Package
5 slide

PRODUCT

P815 (mast cell; mastocytoma)

Summary

Mouse P815 cells were cultured in DMEM and harvested at the log phase of growth. The P815 cells were fixed in 4% paraformaldehyde and arrayed on a 12-well (5 mm) adhesive coated slide, with each well's surface specifically treated to enhance cellular attachment and to minimize background staining. Each well was loaded with approximately 7.5×10^3 P815 cells. A follow-up heat dehydration process was done to ensure the attachment of P815 cells and the stability of cellular proteins.

Applications

Application Note

P815 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

Properties

Storage Stable for up to 6 months at 4°C under desiccate conditions.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).