

Human Heart Ventricle, left (Arrhythmia, infarct) tissue lysate

Cat. No. GTX25576

Species	Human
---------	-------

Package	200 µg
---------	--------

PRODUCT

Summary

Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.

Tissue type

Cardiovascular Disease

Properties

Form	Liquid
------	--------

Buffer	HEPES, MgCl ₂ , KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---

Concentration	5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).