

PTBP1 antibody, N-term

Cat. No. GTX25642

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.01-0.3µg/ml
IHC-P	2-4µg/ml

Note : Human Kidney shows nuclear membrane staining in all the convoluted tubules.

Not tested in other applications.

Calculated MW	57 kDa. (Note)
Product Note	This antibody is expected to recognise all 4 reported isoforms (represented by NP_002810.1; NP_114367.1; NP_114368.1; NP_787041.1).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence DGIVPDIAVGTKR-C, from the N Terminus of the protein sequence according to NP_002810.1; NP_114367.1; NP_114368.1; NP_787041.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX25642 WB Image

WB analysis of mouse ovary (A) and rat ovary (B) lysate using GTX25642 PTBP1 antibody, N-term.

Dilution : 0.1-0.3µg/ml

Loading : 35µg protein in RIPA buffer




GTX25642 IHC-P Image

IHC-P analysis of human kidney using GTX25642 PTBP1 antibody, N-term.

Antigen retrieval : Tris/EDTA buffer pH 9

Dilution : 2µg/ml



GTX25642 WB Image

WB analysis of human ovary lysate using GTX25642 PTBP1 antibody, N-term.

Dilution : 0.1µg/ml

Loading : 35µg protein in RIPA buffer



For full product information, images and publications, please visit our [website](https://www.genetex.com).