

BCMA antibody

Cat. No. GTX25972

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Suggested dilution	Recommended dilution
WB	2-5 µg/ml

Not tested in other applications.

Calculated MW 20 kDa. ([Note](#))

Properties

Form Liquid

Buffer PBS

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Immunogen A synthetic peptide corresponding to the carboxy terminus of human BCMA.

Purification Purified by affinity chromatography

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).