

SNX5 antibody, C-term

Cat. No. GTX25983

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-3µg/ml
IHC-P	3-6µg/ml

Note : Human Breast shows cytoplasm staining of ductal epithelial cells. Cytoplasm staining is also observed in PCT of Human Kidney and in trophoblasts of Human Placenta .

Not tested in other applications.

Calculated MW	47 kDa. (Note)
Product Note	This antibody is expected to recognize both reported isoforms (NP_055241.1; NP_001269383.1). Reported variants represent identical protein: NP_055241.1, NP_689413.1

Properties	
Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence SLLQSCIDLFKNN, from the C Terminus of the protein sequence according to NP_055241.1; NP_001269383.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx25983 IHC-P Image

IHC-P analysis of human breast using GTx25983 SNX5 antibody, C-term.
Antigen retrieval : citrate buffer pH 6
Dilution : 4µg/ml



GTx25983 WB Image

WB analysis of Jurkat lysate using GTx25983 SNX5 antibody, C-term.
Dilution : 1µg/ml
Loading : 35µg protein in RIPA buffer



For full product information, images and publications, please visit our [website](#).