

PHAP antibody

Cat. No. GTX25987

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA, IHC
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/ml
IHC-P	Assay dependent
ELISA	Assay dependent
IHC	10µg/ml

Not tested in other applications.

Calculated MW 29 kDa. ([Note](#))

Product Note This antibody recognizes human, mouse and rat PHAP I (32kDa) and PHAP III (35kDa).

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to amino acids near the carboxy terminus of human PHAP (accession no. P39687).
Purification	Purified by DEAE chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTX25987 IHC-P Image

IHC-P analysis of human small intestine tissue using GTX25987 PHAP antibody.

Dilution : 20 µg/ml


GTX25987 WB Image

Western blot analysis of PHAP expression in Raji cell lysate with PHAP antibody (GTX25987) at 1 µg/ml. The two bands below PHAP I might be differently spliced isoforms of PHAP.


GTX25987 IHC Image

Immunohistochemical staining of PHAP in human small intestine with PHAP antibody (GTX25987) at 10 µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).