

ANP32A antibody

Cat. No. GTX25991

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/ml
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	29 kDa. (Note)
Product Note	This antibody does not cross-react with PHAP I2a or PHAP III.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to amino acids at the carboxy terminus of human PHAP (accession no. P39687).
Purification	Purified by DEAE chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX25991 IHC-P Image

Immunohistochemical staining of PHAP I in mouse small intestine with PHAP I antibody (GTX25991) at 2µg/ml.



GTX25991 IHC-P Image

IHC-P analysis of Mouse Small Intestine cells using GTX25991 ANP32A antibody.
Dilution : 20 µg/ml



GTX25991 WB Image

Western blot analysis of PHAP I expression in human Raji cell (A), mouse (B) and rat (C) testis tissue lysates with PHAP I antibody (GTX25991) at 1µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).