

Goat Anti-Rabbit IgG (Fc fragment) antibody, F(ab')<sub>2</sub> fragment (TxRd)

Cat. No. GTX26019

Host	Goat
Clonality	Polyclonal
Isotype	IgG F(ab') <sub>2</sub>
Applications	ICC/IF, FCM, Dot, ELISA
Reactivity	Rabbit

Package

500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:1000-1:5000
FCM	1:500-1:2500
Dot	Assay dependent
ELISA	1:10000-1:50000

Not tested in other applications.

## Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl, 1% BSA
Preservative	0.01% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit IgG Fc fragment
Purification	IgG fraction This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorptions to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
Conjugation	Texas Red (TxRd) <a href="#">Wavelength</a> Ratio : 4.0 molecules TxRd per Goat IgG F(ab') <sub>2</sub> molecule.

For full product information, images and publications, please visit our [website](#).


Date 2025 / 12 / 27 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTx26019 Dot Image

Dot blot analysis of rabbit IgG Fc fragment using GTx26019 Goat Anti-Rabbit IgG (Fc fragment) antibody, F(ab')<sub>2</sub> fragment (TxRd).

Lane 1 : 100 ng

Lane 2 : 33.3 ng

Lane 3 : 11.1 ng

Lane 4 : 3.7 ng

Lane 5 : 1.23 ng

Dilution : 1 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).