

Rb antibody

Cat. No. GTX26075

Host	Rabbit	
Clonality	Polyclonal	
Isotype	IgG	
Applications	WB	
Reactivity	Human, Mouse, Rat	

Package 1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution		Recommended dilution
WB		Assay dependent
Not tested in other app	olications.	
Calculated MW	106 kDa. (<u>Note</u>)	

Properties		
Form	Liquid	
Buffer	PBS, 0.2% Gelatin	
Preservative	0.1% Sodium azide	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)	
Immunogen	Peptide mapping to the carboxy terminal domain of human Rb.	
Purification	Purified by affinity chromatography	
Conjugation	Unconjugated	
	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or	



Note

For full product information, images and publications, please visit our <u>website</u>.

human consumption.

determine the structure or sequence of the product.

Date 2025 / 11 / 03 Page 1 of 1

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to