

CIKS peptide

Cat. No. GTX26138

Applications	Neutralizing /Inhibition
Species	Human

Package
50 µg

Applications

Application Note

WB: GTX26138 will completely block antibody binding in Western blot by incubating the peptide with an equal concentration of antibody for 30 minutes at 37 °C. Optimal dilutions/concentrations should be determined by the end user.

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store the peptide at -20°C, stable for one year.
Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).