

NAK / TBK1 peptide

Cat. No. GTX26236

Applications	Neutralizing /Inhibition
Species	Human

Package
50 µg

Applications

Application Note

The peptide is used for blocking the antibody activity of anti-NAK (cat. no 8121). It usually blocks the antibody activity completely in Western blot by incubating the peptide with an equal volume of antibody for 30 min at 37°C.

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store the peptide at -20°C, stable for one year.
Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).