

## Human Brain: Parietal Lobe (Parkinson's Disease) tissue lysate

Cat. No. GTX26370

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

|         |        |
|---------|--------|
| Package | 200 µg |
|---------|--------|

## PRODUCT

|         |                                                                                                                                                                                                  |
|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Summary | Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail. |
|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|             |                      |
|-------------|----------------------|
| Tissue type | Neurological Disease |
|-------------|----------------------|

## Properties

|      |        |
|------|--------|
| Form | Liquid |
|------|--------|

|        |                                                                                                                        |
|--------|------------------------------------------------------------------------------------------------------------------------|
| Buffer | HEPES, MgCl <sub>2</sub> , KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors |
|--------|------------------------------------------------------------------------------------------------------------------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                   |
|---------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|---------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                                                                          |
|---------------|--------------------------------------------------------------------------|
| Concentration | 5 mg/ml (Please refer to the vial label for the specific concentration.) |
|---------------|--------------------------------------------------------------------------|

|      |                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Note | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------|

|      |                                                                                                                                                                             |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Note | Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|



For full product information, images and publications, please visit our [website](#).

Date 2026 / 01 / 13 Page 1 of 1