

CDK1 + CDK2 antibody

Cat. No. GTX26434

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

ELISA: Use at an assay dependent dilution. WB: Use at a concentration of 2 µg/ml. Detects a band of approximately 33 kDa in HeLa nuclear extract. Optimal dilutions/concentrations should be determined by the end user.

Product Note Recognises both cdk1 and cdk2 - useful as a general cdk antibody since it reacts with both cdk1 and cdk2 and is very sensitive giving an intense signal.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.7 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human cdk2
Purification	Protein A purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).