

Rabbit Anti-Sheep IgG antibody (AP)

Cat. No. GTX26748

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	WB, Dot, ELISA, IHC
Reactivity	Sheep

Package
1 mg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
Dot	Assay dependent
ELISA	1:5000
IHC	1:200-1:1000

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	50mM Tris-HCl, 150mM NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂ , 1% BSA, 50% Glycerol
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Sheep IgG whole molecule
Purification	IgG fraction This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorptions to remove any unwanted reactivities.
Conjugation	Alkaline Phosphatase (AP)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX26748 Dot Image

Dot blot analysis of sheep IgG using GTX26748 Rabbit Anti-Sheep IgG antibody (AP).



For full product information, images and publications, please visit our [website](#).