

## Rabbit Anti-Guinea Pig IgG antibody (FITC)

**Cat. No. GTX26767**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	ICC/IF, FCM, ELISA
<b>Reactivity</b>	Guinea pig

**Package**  
1 mg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:1000-1:5000
FCM	1:500-1:2500
ELISA	1:10000-1:50000

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	10mM Sodium Phosphate, 150mM NaCl, 1% BSA
<b>Preservative</b>	0.01% Thimerosal
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
<b>Concentration</b>	2.2 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Guinea Pig IgG whole molecule
<b>Purification</b>	Purified by antigen-affinity chromatography. From serum
<b>Conjugation</b>	Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a> Ratio : 2.0 molecules FITC per Rabbit IgG molecule.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).