

Goat Anti-Golden Syrian Hamster IgG antibody (HRP)

Cat. No. GTX26892

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Application	WB, ELISA, IHC
Reactivity	Golden Syrian Hamster

Reference (1)
Package
1 mg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:5000
ELISA	1:7000-1:35000
IHC	1:500-1:2500

Not tested in other applications.

Product Note Reduced reactivity with Armenian Hamster IgG.

PROPERTIES

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl, 1% BSA
Preservative	0.01% Gentamicin Sulfate
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Golden Syrian Hamster IgG whole molecule
Purification	IgG fraction This product was prepared from monospecific antiserum by immunoaffinity chromatography using Hamster IgG coupled to agarose beads followed by solid phase adsorptions to remove any unwanted reactivities.
Conjugation	Horseradish peroxidase(HRP)




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX26892 ELISA Image

ELISA analysis of Golden Syrian Hamster IgG using serially diluted GTX26892 Goat Anti-Golden Syrian Hamster IgG antibody (HRP).

Coating : 1 μ g



For full product information, images and publications, please visit our [website](https://www.genetex.com).