

Goat Anti-Human IgG antibody, pre-adsorbed (Cy2)

Cat. No. GTX26956

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P, IHC-Fr, FCM
Reactivity	Human

Package 100 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100 - 1:1,000
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	1:200 - 1:400

Not tested in other applications.

Product Note

This antibody has been pre-adsorbed against Mouse, Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Rabbit, Rat and Sheep serum proteins to ensure minimal cross-reactivity. Cross-reactivity with other species, however, is possible.

Properties	
Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl, 1% BSA
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG whole molecule
Purification	Immunogen affinity purified
Conjugation	Cyanine2 (Cy2) Wavelength



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 2



For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 2 of 2