

# Goat Anti-Human IgG antibody, pre-adsorbed (Cy2)

**Cat. No. GTX26956**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Application</b>	ICC/IF, IHC-P, IHC-Fr, FACS
<b>Reactivity</b>	Human

**Package**  
100 µg

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100 - 1:1,000
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FACS	1:200 - 1:400

Not tested in other applications.

### Product Note

This antibody has been pre-adsorbed against Mouse, Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Rabbit, Rat and Sheep serum proteins to ensure minimal cross-reactivity. Cross-reactivity with other species, however, is possible.

## PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	0.02M Potassium Phosphate, 0.15M NaCl, 1% BSA
<b>Preservative</b>	0.01% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Human IgG whole molecule
<b>Purification</b>	Immunogen affinity purified
<b>Conjugation</b>	Cyanine2 (Cy2)



For full product information, images and publications, please visit our [website](#).

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).