

NIH-3T3 whole cell lysate

Cat. No. GTX27179

Species Mouse

Package
100 µg

PRODUCT

Summary

NIH-3T3 lysate was prepared by homogenization in modified RIPA buffer. Cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The lysate was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Properties

Form Liquid

Buffer 1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).