

MMP2 antibody

Cat. No. GTX27298

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Rat

Package
50 µg

Applications

Application Note

For IP: Use at a dilution of 1:500. For WB: Use at a dilution of 1:2000. Optimal dilutions/concentrations should be determined by the researcher

Calculated MW	74 kDa. (Note)
Product Note	MMP2 - No cross-reactivity with other MMP family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E. coli-expressed active rat 72-kDa type IV collagenase (catalytic domain)
Purification	Protein A affinity purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx27298 WB Image


For full product information, images and publications, please visit our [website](https://www.genetex.com).