

Protein A (AP)

Cat. No. GTX27345-pro

Species Staphylococcus aureus

Package
250 µg

APPLICATION

Product Note No reaction was observed against anti-Protein G.

PROPERTIES

Form Liquid

Buffer 0.05M Tris Chloride (pH8.0), 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 10mg/ml BSA (IgG, Protease free), 50% (v/v) Glycerol

Preservative 0.01% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Expression System Native protein from Staphylococcus aureus

Conjugation Alkaline Phosphatase (AP)

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).