

## Protein A (AP)

Cat. No. GTX27345-pro

Species	Staphylococcus aureus
---------	-----------------------

Package	250 µg
---------	--------

## Applications

Product Note	No reaction was observed against anti-Protein G.
--------------	--

## Properties

Form	Liquid
------	--------

Buffer	0.05M Tris Chloride (pH8.0), 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 10mg/ml BSA (IgG, Protease free), 50% (v/v) Glycerol
--------	--

Preservative	0.01% Sodium azide
--------------	--------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
---------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Expression System	Native protein from Staphylococcus aureus
-------------------	---

Conjugation	Alkaline Phosphatase (AP)
-------------	---------------------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---



For full product information, images and publications, please visit our [website](#).