

Acinus antibody

Cat. No. GTX27354

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-10µg/ml
ICC/IF	5-10µg/ml
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 152 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide corresponding to aa 775- 789 of human AcinusL, aa 17-31 of human Acinus S and aa 48-62 of human AcinusS' which differ from those of mouse Acinus by one amino acid.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX27354 WB Image

Western blot analysis of Acinus in K562 whole cell lysate with Acinus antibody (IN) (GTX27354) at 1µg/ml.




GTX27354 ICC/IF Image

ICC/IF analysis of K562 cells using GTX27354 Acinus antibody.

Green : Primary antibody

Blue : DAPI

Dilution : 20 µg/ml



GTX27354 IHC-P Image

Immunocytochemical staining of K562 cells using Acinus antibody (GTX27354) at 5µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).