

Interferon gamma antibody

Cat. No. GTX27361

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	ELISPOT
Reactivity	Rat

Package
100 µg

APPLICATION

Application Note

The use of this product can be used in ELISPOT at an assay dependent dilution, but it has not yet been reported. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 18 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein - rat IFNG protein
Purification	Protein A affinity purified From polyclonal serum
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).