

# Bovine Serum (Calf) - Sterile

**Cat. No. GTX27479**

<b>Application</b>	WB, ICC/IF, ELISA, Blocking, IHC
<b>Species</b>	Bovine

**Package**  
25 ml

## PRODUCT

### Summary

Normal Bovine Calf Serum are obtained from non-demolished blood that is collected from healthy, fasted donors. The donors are not medicated in any way and are maintained on an antibiotic free diet. Sterilization is accomplished by Millipore filtration. These serums provide excellent growth promoting components for tissue culture and microbiological organisms.

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
ELISA	Assay dependent
Blocking	Assay dependent
IHC	Assay dependent

Not tested in other applications.

## PROPERTIES

**Form** Liquid

**Storage** Keep as concentrated solution. Store at -20°C or below. Avoid multiple freeze-thaw cycles. Use standard sterile procedures when handling this product to maintain sterility.

**Purification** Sterilization: This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).