

Goat Serum - Sterile

Cat. No. GTX27481

Applications	ICC/IF
Species	Goat

References (25)
 Package
 50 ml

PRODUCT

Summary Goat Serum

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

ICC/IF	Assay dependent
--------	-----------------

Not tested in other applications.

Properties

Form Liquid

Buffer Serum


Preservative No preservatives

Storage Keep as concentrated solution. Store at -20°C or below. Avoid multiple freeze-thaw cycles. Use standard sterile procedures when handling this product to maintain sterility.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX27481 Image

The data was published in the journal Cell Death Dis in 2019. [PMID: 31127082](https://pubmed.ncbi.nlm.nih.gov/31127082/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).



GTX27481 Image

The data was published in the journal Cell Death Dis in 2019. [PMID: 31127082](https://pubmed.ncbi.nlm.nih.gov/31127082/)



For full product information, images and publications, please visit our [website](#).