

Bovine Collagen IV protein

Cat. No. GTX27529-pro
Species Bovine

Package
500 µg

Applications

Observed MW (kDa) 48 (Collagen IV alpha 3), 58 (Collagen IV alpha 2), 65 (Collagen IV alpha 1) kDa.

Properties

Form Liquid

Buffer 0.5M Acetic acid

Preservative 0.01% Sodium azide

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.


Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Expression System Bovine Placenta

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX27529-pro Image

SDS-PAGE analysis of GTX27529-pro Bovine Collagen IV protein.

Lane 1 : Bovine Collagen Type IV Reduced (5.0µg)

Lane 2 : Prestained Molecular Weight Marker

Lane 3 : Bovine Collagen Type IV Non-Reduced (5.0µg)



For full product information, images and publications, please visit our [website](https://www.genetex.com).