

## Bovine Collagen IV protein

Cat. No. GTX27529-pro

Species	Bovine
---------	--------

Package	500 µg
---------	--------

## Applications

**Observed MW (kDa)** 48 (Collagen IV alpha 3), 58 (Collagen IV alpha 2), 65 (Collagen IV alpha 1) kDa.

## Properties

Form	Liquid
------	--------

Buffer	0.5M Acetic acid
--------	------------------

Preservative	0.01% Sodium azide
--------------	--------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---


Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Expression System	Bovine Placenta
-------------------	-----------------

Conjugation	Unconjugated
-------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

## DATA IMAGES



## GTX27529-pro Image

SDS-PAGE analysis of GTX27529-pro Bovine Collagen IV protein.

Lane 1 : Bovine Collagen Type IV Reduced (5.0µg)

Lane 2 : Prestained Molecular Weight Marker

Lane 3 : Bovine Collagen Type IV Non-Reduced (5.0µg)



For full product information, images and publications, please visit our [website](#).