

## Bovine Collagen V protein

Cat. No. GTX27530-pro

Species	Bovine
---------	--------

Package	500 µg
---------	--------

## Properties

Form	Liquid
------	--------

Buffer	0.2M Acetic acid (pH2.5)
--------	--------------------------

Preservative	0.01% Sodium azide
--------------	--------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Expression System	Bovine Placenta
-------------------	-----------------

Conjugation	Unconjugated
-------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---



For full product information, images and publications, please visit our [website](#).