

# Bovine Collagen V protein

**Cat. No. GTX27530-pro**

<b>Species</b>	Bovine
----------------	--------

<b>Package</b> 500 µg
--------------------------

## Properties

<b>Form</b>	Liquid
-------------	--------

<b>Buffer</b>	0.2M Acetic acid (pH2.5)
---------------	--------------------------

<b>Preservative</b>	0.01% Sodium azide
---------------------	--------------------

<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	---

<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	--

<b>Expression System</b>	Bovine Placenta
--------------------------	-----------------

<b>Conjugation</b>	Unconjugated
--------------------	--------------

<b>Note</b>	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---



For full product information, images and publications, please visit our [website](#).