

Human Collagen I protein


Cat. No. GTX27533-pro**Species** Human**Package**

500 µg

Properties

Form Liquid**Buffer** 0.5M Sodium Acetate (pH4.5)**Preservative** 0.01% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Expression System** Human Placenta**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES

**GTX27533-pro Image**

Western blot of Human Collagen Type I (GTX27533-pro) using Collagen Type I antibody.

Load: 50 ng.

Primary antibody: Collagen Type I antibody, 1:1,000 overnight at 4°C.

Predicted/Observed size: 139 & 130 kDa, 139 & 130 kDa for Collagen Type I.

Other Band(s): Collagen Type I splice variants and isoforms.

For full product information, images and publications, please visit our [website](#).