

Human Collagen V protein

Cat. No. GTX27537-pro

Species	Human
----------------	-------

Package 500 µg

Properties

Form	Liquid
Buffer	0.5M Acetic acid
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Expression System	Human Placenta
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).