

IκB alpha peptide

Cat. No. GTX27546

Applications	Neutralizing/Inhibition
Species	Human

Package
50 µg

Applications

Application Note

Intended for use as a control peptide when used with anti-IκB alpha (C-terminal specific ab GTX27545) to block specific interaction of anti-IκB alpha (C-terminal specific) with the IκB alpha subunit. Control peptide should be used at 1.0 mg per 1.0 ml of antiserum per assay.

Properties

Form	Liquid
Buffer	0.15M Sodium Chloride, 0.02M Potassium Phosphate pH7.2
Preservative	0.01% Sodium azide
Storage	Keep as concentrated solution. Store at 4°C short term. For extended storage aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).