

IKB beta antibody

Cat. No. GTX27547

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Mouse, Rat

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ELISA	1:5000-1:25000

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Filtered-steriled serum
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	IkBb peptide corresponding to a region near the C-terminus of the human protein conjugated to KLH.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX27547 WB Image

WB analysis of various samples using GTX27547 IKB beta antibody.

Lane 1 : Protein ladder

Lane 2 : A549 whole cell lysate

Lane 3 : A431 whole cell lysate

Lane 4 : Molt-4 whole cell lysate

Loading : 10 _x0007_μL

Dilution : 1:500



For full product information, images and publications, please visit our [website](#).