

IκB beta peptide

Cat. No. GTX27548

Applications	Neutralizing /Inhibition
Species	Human

Package
50 µg

Applications

Application Note

For BL: Intended for use as a control peptide when used with anti-IκB beta (C-terminal specific GTX27548) to block specific interaction of anti-IκB beta (C-terminal specific) with the IκB beta subunit. Control peptide should be used at 1.0 mg per 1.0 ml of antiserum in per assay. Optimal dilutions/concentrations should be determined by the researcher.

Properties

Form	Liquid
Buffer	Batch dependent
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purity	> 95%
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 1