

# Human Parotid (Tumor) tissue lysate, Membrane Fraction

**Cat. No. GTX27677**

**Species** Human

**Package**  
100 µg

## PRODUCT

### Summary

Human parotid tissue membrane protein lysate was prepared by isolating the membrane protein from whole tissue homogenates using a proprietary technique. The human parotid was frozen in liquid nitrogen immediately after excision and then stored at -70°C. For quality control purposes, the isolated parotid tissue membrane protein pattern on SDS-PAGE gel is shown to be consistent for each lot by visualization with coomassie blue staining.

### Tissue type

Tumor

## Properties

**Form** Liquid

**Buffer** HEPES, MgCl<sub>2</sub>, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

**Preservative** No preservatives

### Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

### Concentration

2 mg/ml (Please refer to the vial label for the specific concentration.)

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).