

Human Pancreas (Diabetes) tissue lysate, Membrane Fraction

Cat. No. GTX27711

Species	Human
---------	-------

Package	100 µg
---------	--------

PRODUCT

Summary	Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.
---------	--

Tissue type	Diabetes
-------------	----------

Properties

Form	Liquid
------	--------

Buffer	HEPES, MgCl ₂ , KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---

Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 16 Page 1 of 1