

# TNF alpha antibody

**Cat. No. GTX27742**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, FCM, Immunoassay
<b>Reactivity</b>	Human, Monkey

**Package**  
250 µl

## Applications

### Application Note

For FACS, IA, and WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

<b>Calculated MW</b>	26 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	The antibody can be used for neutralization of human and monkey TNF-alpha in bio-assays.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 0.1% BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Purification</b>	IgG fraction
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).