

TNF alpha antibody

Cat. No. GTX27742

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	WB, FACS, Immunoassay
Reactivity	Human, Monkey

Package
250 µl

APPLICATION

Application Note

For FACS, IA, and WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	26 kDa. (Note)
Product Note	The antibody can be used for neutralization of human and monkey TNF-alpha in bio-assays.

PROPERTIES

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).