

alpha Tubulin antibody [TU-16]

Cat. No. GTX27749

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	WB, ICC/IF, IHC-P, IP, ELISA
Reactivity	Human, Mouse, Rat, Bovine, Hamster, Pig

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties


Form	Liquid
Buffer	TBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Porcine brain microtubule protein MTP-1.
Purification	Purified by precipitation and chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX27749 IP Image

IP analysis of HeLa and Raji cell lysate using GTX27749 alpha Tubulin antibody [TU-16]. The alpha tubulin signal was detected by GTX27750 alpha Tubulin antibody [TU-01]

L : Cell lysate

IPr : Immunoprecipitate lysate (reducing conditions)

IPn : Immunoprecipitate lysate (non-reducing conditions)



For full product information, images and publications, please visit our [website](https://www.genetex.com).