

NF-M antibody [NF-09]

Cat. No. GTX27794

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat, Bovine, Pig

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Product Note

This antibody reacts with both phosphorylated and non-phosphorylated form of medium neurofilament protein (160 kDa) of various species.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pellet of porcine brain cold stable proteins after depolymerization of microtubules.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX27794 WB Image**

WB analysis of pig brain tissue lysate using GTX27794 NF-M antibody [NF-09].

**GTX27794 ICC/IF Image**

ICC/IF analysis of Neuro2a cells using GTX27794 NF-M antibody [NF-09].

Red : Primary antibody

Blue : DAPI



For full product information, images and publications, please visit our [website](#).