

CTGF peptide

Cat. No. GTX27861

Applications	Neutralizing/Inhibition
Species	Mouse

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing/Inhibition	Assay dependent
-------------------------	-----------------

Not tested in other applications.

Properties

Storage	Keep as concentrated solution. Store at 4°C short term. For extended storage aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	---

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	--

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
-------------	--



For full product information, images and publications, please visit our [website](#).