

Human Lymphoma, Hodgkins Disease (Tumor) tissue lysate, Membrane Fraction

Cat. No. GTX27889

Species Human**Package**

100 µg

PRODUCT

Summary

Human lymphoma, Hodgkins Disease tissue membrane protein lysate was prepared by isolating the membrane protein from whole tissue homogenates using a proprietary technique. The human lymphoma, Hodgkins Disease tissue was frozen in liquid nitrogen immediately after excision and then stored at -70°C. For quality control purposes, the isolated lymphoma, Hodgkins Disease tissue membrane protein pattern on SDS-PAGE gel is shown to be consistent for each lot by visualization with coomassie blue staining.

Tissue type

Tumor

Properties

Form

Liquid

Buffer

HEPES, MgCl₂, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

Preservative

No preservatives

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration

2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).