

Human Brain (Tumor) tissue lysate

Cat. No. GTX27925

Species Human**Package**
100 µg

PRODUCT

Summary

Brain tissue lysate was prepared by homogenization in modified RIPA buffer. Tissue and cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Tissue type

Tumor


Properties

Form Liquid**Buffer** 1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 2 mg/ml (Please refer to the vial label for the specific concentration.)**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX27925 Image



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 1