

NFkB p105/p50 antibody

Cat. No. GTX27971

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, IHC
Reactivity	Human, Mouse, Rat

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:400
IP	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 105 kDa. ([Note](#))

Product Note This antibody reacts with NFkB p50 and p105.

Properties

Form	Liquid
Buffer	PBS, 0.2% Gelatin
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide mapping to the carboxy terminus of human NFkB p50.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX27971 WB Image

WB analysis of HeLa (A), HL-60 (B), and NIH/3T3 (C) whole cell lysates using GTX27971 NFkB p105/p50 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).